

[HMK-200 Air Jet Sieve Specifications](#)

- ◇ Nozzle Rotation: 35 rpm
- ◇ Nozzle Type: 2mmx1, 1mmx2
- ◇ Sieving Time: Timed set up to 99m 59s
- ◇ Dimensions: W35xD35xH35cm
- ◇ Power: 220V/110V 50/60Hz 25W
- ◇ Weight: 11.8Kg
- ◇ Sieving Range: 5 to 5,000 micron
- ◇ Sample Size: 0.1 to 2,000 grams

[HMK-200 Air Jet Sieving Machine](#)