

**-Flexible parameter adjustment**

The AS-100 Tapped Density Tester changes easily the key parameter such as tapping speed and drop height in order to meet various international standards.

**-Digital LCD Display, Fast and easy operation**

The LED and membrane control panel brings operator an easy and fast operation. Accurate to 4 decimals. Automatically saving all the test information.

**-Reduced noise**

AS-100 is equipped with a safeguard acoustic cabinet in order to reduce the noise level, while operator can monitor the operation through the front transparent plexiglass door.

**-Easily outputting hard copy**

The integrated micro computer prints the test data including Sample No, Weight, Volume, Tap Number, Tap Density and Operator by simply pressing the button on membrane.

**-Easy to clean**

The outer surface is plastic coated, it is strong and resistant to any powders non-corrosive or corrosive. It makes for the cleaning procedure as well as a longer working life.